




Changzhou Haida Medical Equipment Co.,Ltd.

**Product name: Disposable Mask**

(for personal use in public places **NOT** for hospital use)

Product Name	Disposable mask
Size	175*95mm
Product picture	
Material	<p>1<sup>st</sup> layer: Hydrophobic non-woven fabric 2<sup>nd</sup> layer :Filter melt blown cloth 3<sup>rd</sup> layer: Skin friendly non-woven fabric</p> 
Packing	<p>50 PCS in a box</p> 

Changzhou Haida Medical Equipment Co.,Ltd.

40boxes (2000pcs) in a carton



Gross weight:9.5-10KGS

Carton size: 520\*400\* 380 mm

Label



In English

(sample)



In Chinese

Instruction For use

See attachment- IFU in English

Chinese IFU below

## 一次性使用口罩说明书

【执行标准】Q/320402 HD 001-2020

【产品结构】一次性使用口罩（以下简称口罩）由热轧拒水无纺布、熔喷布、鼻夹、口罩带组成。

【产品规格】成人型17.5cm\*9.5cm

【适用范围】一次性使用口罩适用于覆盖使用者的口、鼻及下颌，用于日常生活中或普通医疗环境中佩戴、阻隔口腔和鼻腔呼出或喷出污染物的一次性口罩。

【使用方法】

1. 取出一一次性使用口罩，将口罩有鼻夹的一边向上，两端的口罩带挂于双耳；
2. 轻按鼻夹使其与鼻梁贴合，然后按着鼻夹，将口罩下端拉至下颌处；
3. 调整口罩佩戴效果，使其覆盖鼻梁至下颌处，以获得最佳的防护效果。

【注意事项】

1. 本产品一次性使用，过期产品请勿使用。
2. 对非织布过敏者慎用。
3. 产品使用后，应按环保部门要求进行处理。
4. 本产品不得代替医用外科口罩或医用防护口罩。

【储存运输】

产品在储运中应保证密封、不破损、不受潮，注意防火、防雨、防酸、防碱、避免强光直射。

【生产日期】见合格证      【有效期】见合格证

生产企业：常州市海达医疗器械有限公司

生产地址：常州市天宁区郑陆镇三河口工业集中区

电话：0519-88679868

传真：0519-88672999

邮编：213115

CMA/CNA test  
report

Available

Ce declaration

Available